

Invacare® | **ProFlex Mattress Overlay**



**1 Year warranty**




**Good value pressure redistributing overlay**

The **ProFlex Overlay** has been designed to be used in conjunction with an existing mattress to provide an effective Pressure Ulcer Management solution. Designed for those at 'At Risk/Increased Vulnerability'\* of developing pressure ulcers, the Proflex overlay is constructed from high quality contoured foam to maximise pressure reduction. The cover is a 2 way stretch waterproof and vapour permeable fabric.



**Yes, you can.®**




**Features and Options**

**Effective Pressure Re-distribution**

The contoured foam affords pressure re-distribution and cradles bony prominences. Made from high quality foam, the **ProFlex Overlay** offers a good life span, and the possibility of the overlay "bottoming out" is reduced. **ProFlex Overlay** sits on top of the existing mattress to provide a comfortable surface for the patient.

**Even Weight Distribution**

**ProFlex Overlay** has cut foam castellations to reduce shear and friction forces. The vapour permeable properties of the cover, mean the micro climate is maintained and there is less risk of sweating and associated tissue viability complications.



**Cover helps to prevent a build up of shear and friction forces**

Cover is multi stretch, water-resistant and vapour-permeable.



**Covers have an anti-slip base**

Ensures safe use of the **Pro Flex Overlay** and helps to prevent slippage.

**Technical data**

<b>ProFlex Mattress Overlay</b>					
	Std sizes: 187.5 x 87.6 x 8.5 cm	Height: 8.5 cm	Weight: 5.5 kg	Max User Weight: 108 kg (17 stone)	Cleaning: 80°

ProFlex Overlay Product Code:  
1539207

Manufactured to European Standards:  
BS 7177 Crib 5

**Warranty 1 year**

Invacare Limited  
Pencoed Technology Park - Pencoed -  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: uk@invacare.com  
Sales Order E-mail:  
ordersuk@invacare.com  
www.invacare.co.uk

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS - County  
Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: Ireland@invacare.com  
www.invacare.ie

© 2015 Invacare International Sàrl  
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

ProFlex Mattress Overlay  
UK - 03/2015 - 1585806

**Colour**



Baltic Blue

\* Essential nursing care is pivotal in pressure ulcer prevention. This product will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.



**Yes, you can.®**